

Conners CPT 3™, Conners CATA™, and Conners K-CPT 2™: Introduction and Application



Amy Patenaude, Ed.S., NCSP
MHS Assessment Consultant ~ School Psychologist




Your Presenter


I put years of studying, sweat
and tears into it,
so yes, I am going
to put NCSP after
my name on
EVERYTHING!



somecards
© 2013

Amy Patenaude, Ed.S., NCSP
MHS Assessment Consultant
Amy.patenaude@mhs.com

 @Amy_Patenaude

 <https://www.linkedin.com/in/amyptenaude>



Disclosure

I work for Multi-Health Systems (MHS),
the publisher of the assessment tools that we will be discussing today.
Rating scales should not be used as the sole basis
for making a diagnosis or educational eligibility decision.



Training Objectives

Best practices in the assessment of ADHD and EF

CPT 3 and CATA and their essential features

CEFI and its essential features

Case study examples

MHS
ASSESSMENTS

CPT3 & CATA: Key Features



MHS
ASSESSMENTS

What is a Continuous Performance Test (CPT)?

- Performance/Task based assessment that measures different areas of attention such as sustained attention, inattentiveness, impulsivity, and vigilance
- Provides **objective** information regarding an individual's attentional difficulties
- Clients presented with repetitive boring task and must maintain their focus over a period of time in order to respond to targets or inhibit response to non-targets

CONNERS CPT3 & CATA

MHS
ASSESSMENTS

Why should we be interested in utilizing the
Conners CPT-3 and Conners CATA?

- Objective measure
- Performance based task: engages students/patients/clients and helps build rapport
- Helps pinpoint type of attention problem
- Increases diagnostic and classification accuracy when paired with other assessment measures

CONNERS
CPT 3
Continuous Performance Test 3rd Edition

CONNERS
CATA
Continuous Attention Test

MHS
ASSESSMENTS

Why should I be interested in the CPT-3 and
CATA?

Classification Statistic	Conners 3-P	Conners 3-P & Conners CATA	Conners 3-P & Conners CPT 3	Conners 3-P Conners CPT 3 & Conners CATA
Overall Correct Classification (%)	83.9	88.4	88.4	93.8
Sensitivity (%)	86.0	91.2	89.5	94.7
Specificity (%)	81.8	85.5	87.3	92.7

CONNERS
CPT 3
Continuous Performance Test 3rd Edition

CONNERS
CATA
Continuous Attention Test

MHS
ASSESSMENTS

Conners Continuous
Performance Test 3rd Edition
(Conners CPT 3™)

CONNERS
CPT 3
Continuous Performance Test 3rd Edition

CONNERS
CATA
Continuous Attention Test

MHS
ASSESSMENTS

Conners CPT-3 Key Features

- Ages 8+; assesses attention related problems
- 14 minutes; in addition to 1 minute practice test
- Non-X paradigm: ignore X and respond to all other targets
- High proportion of targets to non-targets
- Varied time intervals between targets (1, 2, or 4 sec ISI)
- By-Block statistics (6 blocks with 60 trials each)
- Practice Test
- Can be part of a battery of assessments for ADHD and other disorders/neurological problems characterized by attention problems

CONNERS
CPT3
CANTAB
ATTENTION TESTS

MHS
ASSESSMENTS

What does the Conners CPT-3 measure?

- Assesses attention related problems
- Examines four dimensions of attention:
 1. Inattentiveness
 2. Impulsivity
 3. Sustained Attention
 4. Vigilance
- Validity Check
- Response Style Analysis:
 1. Liberal
 2. Conservative
 3. Balanced

CONNERS
CPT3
CANTAB
ATTENTION TESTS

MHS
ASSESSMENTS

CPT 3 Scores

Variable	Description
C	Assesses Response Style
d'	Ability to discriminate between targets (non-X) and non-targets(X)
Omissions	Misued targets (non-X)
Commissions	Incorrect responses to non-targets(X)
Perserations	Random, repetitive, or anticipatory responses (i.e., HRT < 100ms)
HRT Reaction Time (HRT)	Response Speed
HRT SD/Variability	Response Speed Consistency
HRT Block Change	Change in HRT across blocks of trials
HRT ISI Change	Change in HRT across ISIs

CONNERS
CPT3
CANTAB
ATTENTION TESTS

MHS
ASSESSMENTS

CPT-3 Scores

Dimension	Score	Description
Inattentiveness	Detectability (d')	Ability to discriminate between targets (non-X) and non-targets (X)
	Omissions	Mixed targets
	Commissions	Incorrect responses to non-targets
	Hit Reaction Time (HRT)	Response speed
Impulsivity	HRT Standard Deviation (SD)	Response speed consistency
	HRT	Response speed
	Commissions	Incorrect responses to non-targets
Sustained Attention	Perseverations	Random or anticipatory responses (i.e., HRT < 100ms)
	HRT Block Change	Change in response speed across blocks of trials
	Omissions by block	Mixed targets by block
Vigilance	Commissions by block	Incorrect responses to non-targets by block
	HRT Inter-Stimulus Interval (ISI) Change	Change in response speed at various ISIs
	Omissions by ISI	Mixed targets by ISI
	Commissions by ISI	Incorrect responses to non-targets by ISI

CONNERS CPT3
Continued Attention

CONNERS CATA
Continued Attention

MHS
ASSESSMENTS

Administration

Hardware and Software Requirements

- Intel Core i3 or equivalent performance (recommended)
- 2 GB Ram
- Windows XP or higher
- 1 available USB port
- 12" monitor or larger with minimum resolution of 1024 x768 pixels
- Wired mouse or keyboard

CONNERS CPT3
Continued Attention

CONNERS CATA
Continued Attention

MHS
ASSESSMENTS

New Auditory Test of Attention

Conners Continuous
Auditory Test of Attention
(Conners CATA™)

CONNERS CPT3
Continued Attention

CONNERS CATA
Continued Attention

MHS
ASSESSMENTS

CATA Key Features

- Ages 8+
- Assesses auditory attention and attention problems
- Can be used on its own or as a compliment to the CPT-3 in an assessment battery
- 14 minutes, 200 scored trials, divided into 4 blocks
- Consists of two basic sounds: a low tone and a high tone
- On 80 percent of the trials, the low tone is played first followed by a high tone (warned trial). High tones on warned trials are the targets (AX paradigm)

CONNERS
CPT-3
CATA

MHS
ASSESSMENTS

CATA-Key Features

- On remaining 20 percent of the trials, a high tone is played alone without the low tone (unwarned trial). High tones on unwarned trials are non-targets.
- On most warned trials, the two tones are played sequentially in same ear (non-switch trial)
- On some warned trials, the two tones are played in opposite ears (switch trials)

CONNERS
CPT-3
CATA

MHS
ASSESSMENTS

What does the Conners CATA measure?

- Assesses *auditory processing and attention-related* problems in individuals aged 8 years and older
- Examines three dimensions of attention:
 1. Inattentiveness
 2. Impulsivity
 3. Sustained Attention
- Examines two dimensions of auditory processing:
 1. Auditory Laterality
 2. Auditory Mobility
- Validity Check
- Response Style Analysis:
 1. Liberal
 2. Conservative
 3. Balanced

CONNERS
CPT-3
CATA

MHS
ASSESSMENTS

Conners CATA Scores

Variable	Description
C	Assesses Response Style
d'	Ability to discriminate targets (warned high tone) from non-targets (unwarned high tone)
Omissions	Missed targets (warned high tone)
Commissions	Responded to non-targets (unwarned high tone)
Persistent Commissions	Responded to low sound/Responded before the high sound
HRT	Hit React Time
HRT SD	Response Speed Consistency
HRT Block Change	Change in HRT across blocks
Laterality	HRT & Hits % Left vs. Right Ear (Preference for left vs. right targets)
Mobility	HRT Switch vs. Switch vs. Non Switch Trials (Ability to switch attention from one ear to another)

CONNERS
CPT3
CATA

MHS
ASSESSMENTS

CATA Scores

Dimension	Score	Description
Inattentiveness	Detectability (d')	Ability to discriminate targets (warned high tone) from non-targets (unwarned high tone)
	Omissions	Missed targets
	Commissions	Incorrect responses to non-targets
	Hit Reaction Time (HRT) HRT Standard Deviation (SD)	Response speed Response speed consistency
Impulsivity	HRT	Response speed
	Commissions	Incorrect responses to non-targets
	Persistent Commissions	Incorrect responses before targets
Sustained Attention	HRT Block Change	Change in response speed across blocks of trials
	Omissions by block	Missed targets by block
	Commissions by block	Incorrect responses to non-targets by block
Auditory Laterality	HRT & Hits Left vs. Right Ear	Preference for left vs. right targets
Auditory Mobility	HRT & Hits on Switch vs. Non Switch Trials	Ability to switch attention from one ear to the other

CONNERS
CPT3
CATA

MHS
ASSESSMENTS

Conners CPT 3 & Conners CATA

Standardization & Basic Psychometrics

CONNERS
CPT3
CATA

MHS
ASSESSMENTS

Sample Descriptions

Conners CPT-3

- Normative Sample:
 - N = 1,400 (700 male, 700 female)
 - Spread across the ages
 - 2010 Census Match: Race, Region, (Parental) Education Level
- ADHD Sample:
 - 259 children, 97 adults
 - 62% male
 - 60% medicated

Conners CATA

- Normative Sample:
 - N = 1,080 (540 male, 540 female)
 - Spread across the ages
 - 2010 Census Match: Race, Region, (Parental) Education Level
- ADHD Sample:
 - 193 cases
 - 64% male
 - 63% children (age 8-17)



Conners CPT 3 Split-half Reliability

Variable Type	Measure	Normative		Clinical	
		Children N = 775-800	Adults N = 591-600	Children N = 114-149	Adults N = 136-145
Detectability	d'	.95	.92	.95	.94
	Omissions	.94	.96	.97	.95
Error Type	Commissions	.94	.91	.92	.95
	Perseverations	.90	.73	.95	.90
Reaction Time Statistics	HRT	.99	.99	.98	.99
	HRT SD	.96	.95	.97	.97
	Variability	.80	.73	.85	.79
	Block Change	.90	.91	.80	.91
	ICI Change	.90	.93	.91	.93
Response Style	C	.87	.83	.89	.92



Conners CPT 3 Test-retest Reliability

Variable Type	Measure	Corr*
Detectability	d'	.74**
	Omissions	.83**
Error Type	Commissions	.85**
	Perseverations	.88**
Reaction Time Statistics	HRT	.89**
	HRT SD	.68**
	Variability	.56**
	HRT Block Change	.12
Response Style	HRT ICI Change	.66**
	C	.63**

*Range restriction corrections applied



CATA Split-half Reliability

Variable Type	Measure	Normative		Clinical	
		Children N = 565- 600	Adults N = 480- 480	Children N = 109-122	Adults N = 66-71
Detectability	d'	.97	.98	.96	.93
	Omissions	.93	.94	.98	.97
Error Type	Commissions	.99	.99	.93	.88
	Perseverations	.99	.99	.99	.99
Reaction Time Statistics	HRT	.91	.93	.98	.99
	HRT SD	.86	.90	.81	.95
	Block Change	.96	.95	.90	.92
Response Style	C	.90	.93	.91	.90



CATA Test-retest Reliability

Variable Type	Measure	Corr*
Detectability	d'	.74**
Error Type	Omissions	.65**
	Commissions	.72**
Reaction Time Statistics	Perseverations	.95**
	HRT	.56**
	HRT SD	.63**
Response Style	HRT Block Change	.12
	C	.14

*Range restriction corrections applied



Conners CPT 3 - Group Differences

ADHD vs General Population

Measure		ADHD	Matched Gen. Pop.	Cohen's d	p
		N = 241-246	N = 345-346		
d'	M	-1.9	-2.3	0.43	< .001
	SD	1.0	1.0		
Omissions	M	4.6	3.1	0.25	.001
	SD	6.2	6.2		
Commissions	M	50.1	43.9	0.35	< .001
	SD	20.3	20.3		
Perseverations	M	1.2	0.5	0.38	< .001
	SD	1.7	1.7		
HRT	M	418.7	416.8	0.10	.186
	SD	50.9	60.9		
HRT SD	M	0.305	0.288	0.49	< .001
	SD	0.095	0.095		
Variability	M	0.085	0.069	0.42	< .001
	SD	0.038	0.038		
HRT Block Change	M	0.008	0.003	0.21	.002
	SD	0.023	0.023		
HRT ISI Change	M	0.068	0.059	0.22	.004
	SD	0.052	0.052		



Conners CATA - Group Differences

ADHD vs General Population

Score		ADHD		Matched General Population		F	p	Cohen's d
		N=143-183	M	N=150-203	M			
d'	M	-2.4	-2.2			36.0	< .001	0.57
	SD	1.5	1.3					
Omissions	M	6.6	6.4			1.2	.274	0.10
	SD	15.1	12.1					
Commissions	M	19.0	9.2			43.9	< .001	0.63
	SD	17.9	16.5					
Perseverative Commissions	M	7.7	3.8			7.6	.006	0.26
	SD	15.5	14.3					
HRT	M	293.6	288.1			2.9	.089	0.16
	SD	282.7	272.6					
HRT SD	M	0.360	0.303			13.9	< .001	0.35
	SD	0.140	0.128					
HRT Block Change	M	0.077	0.097			1.8	.182	0.13
	SD	0.082	0.077					

CONNERS
CPT3
CATA

MHS
ASSESSMENTS

7 Step Interpretation Process

- Step 1: Determine Validity of the Administration
- Step 2: Review Response Style Analysis
- Step 3: Examine the Overview of Scores
- Step 4: Review the Overall Summary and Clinical Likelihood
- Step 5: Examine the Individual Dimensions of Attention
- Step 6: Integrate Results with Multiple Sources
- Step 7: Report Results

CONNERS
CPT3
CATA

MHS
ASSESSMENTS

CASE STUDY: GRANT



CONNERS
CPT3
CATA

MHS
ASSESSMENTS

Case Study

- Grant S.
- 10-year-old boy
- Fell behind in school work
- Often seemed distracted and had problems remembering learned materials
- Some ADHD in family history
- Tested for attention deficits using CPT 3 and CATA



MHS ASSESSMENTS

Assessment Plan

1. Clinician to review all available information
2. Obtain primary and differential diagnosis as well as to establish a general picture of Grant's mental and overall health status
3. Administer the following assessments: Conners CPT 3, Conners –March Developmental Questionnaire (CMDQ), Conners 3rd Edition (Conners 3-Parent, Teacher, Self), Conners CATA
4. Systematic clinical interview



MHS ASSESSMENTS

Case Study

Test	Score	Interpretation
Conners CPT 3	10	Normal
Conners 3-Parent	15	Normal
Conners 3-Teacher	12	Normal
Conners 3-Self	11	Normal
Conners CATA	10	Normal

MHS ASSESSMENTS

Case Study

- Step 1: Validity of Administration

Validity of Administration

The Conners CPT 3 performs a validity check based on the number of hits and omission errors committed, as well as a self-diagnostic check of the accuracy of the timing of each administration. If there is an insufficient number of hits to compute scores, and/or if the omission error rate exceeds 25%, these issues will be noted. Also, the program will issue a warning message noting that the administration was invalid if a timing issue is detected.

There was no indication of any validity issues; the current administration should be considered valid.

CONNERS
CPT 3
CATA

MHS
ASSESSMENTS

Case Study

- Step 2: Response Style Analysis


Response Style Analysis

The variable C represents an individual's natural response style in trials that involve a speed-accuracy trade-off. Based on hits or hit score on this variable, a respondent can be classified as having one of the following three response styles: a **conservative style** (T-score ≥ 40) of responding that emphasizes accuracy over speed; a **liberal style** (T-score ≤ 40) of responding that emphasizes speed over accuracy; or a **balanced style** (T-score = 41-50) of responding that is sensitive to both speed and accuracy. Based on Grant's response, he has a **conservative style of responding that emphasizes accuracy over speed** (T-score = 48, 95% Confidence Interval = 34-66). This response style is often associated with slower reaction times, more commission errors (failure to respond to targets), and fewer commission errors (incorrect responses to non-targets). The influence of Grant's conservative response style on other Conners CPT 3 scores should be taken into consideration throughout the interpretation process.

CONNERS
CPT 3
CATA

MHS
ASSESSMENTS

Step 3: Examine the overall profile



CONNERS
CPT 3
CATA

MHS
ASSESSMENTS

Step 4: Clinical Likelihood Statement

Summary: Relative to the normative sample, Grant was less able to differentiate targets from non-targets, made more omission errors, made more commission errors, made more perseverative errors, responded more slowly, displayed less consistency in response speed and displayed more variability in response speed.

Overall, Grant has a total of 7 atypical T-scores, which is associated with a high likelihood of having a disorder characterized by attention deficits, such as ADHD. Note that other psychological and/or neurological conditions with symptoms of impaired attention can also lead to atypical scores on the Conners CPT-3.

Grant's profile of scores and response pattern indicates that he may have issues related to:

- Inattentiveness (Strong Indication)




Step 5: Examine the Individual Dimensions of Attention






Step 5: Examine the Individual Dimensions of Attention





Step 5: Examine the Individual Dimensions of Attention




MHS is the most widely used of attention tests for the adult population. It was first developed in 1973 by Dr. Robert L. MHS, a clinical psychologist at the University of Minnesota. The test is designed to measure the ability to sustain attention and to detect changes in the level of attention. It is a measure of the ability to sustain attention and to detect changes in the level of attention. It is a measure of the ability to sustain attention and to detect changes in the level of attention.

CONNER'S
CPT3
CATA

MHS
ASSESSMENTS

Step 5: Examine the Individual Dimensions of Attention




MHS is the most widely used of attention tests for the adult population. It was first developed in 1973 by Dr. Robert L. MHS, a clinical psychologist at the University of Minnesota. The test is designed to measure the ability to sustain attention and to detect changes in the level of attention. It is a measure of the ability to sustain attention and to detect changes in the level of attention.

CONNER'S
CPT3
CATA

MHS
ASSESSMENTS

Step 5: Examine the Individual Dimensions of Attention



MHS is the most widely used of attention tests for the adult population. It was first developed in 1973 by Dr. Robert L. MHS, a clinical psychologist at the University of Minnesota. The test is designed to measure the ability to sustain attention and to detect changes in the level of attention. It is a measure of the ability to sustain attention and to detect changes in the level of attention.

CONNER'S
CPT3
CATA

MHS
ASSESSMENTS

Step 6: Integrate Results from Multiple Sources

- **CPT3:** problems with inattentiveness
- **CMDQ:** Grant's Uncle diagnosed with ADHD.
- **Conner 3 (P, T, S):** Results suggest problems with inattention. Impairment items related to schoolwork/grades were endorsed.
- **Conners CATA:** problems with inattentiveness
- **Interview:** difficult to get Grant to conduct homework, careless when following instructions. Described Grant as shy and anxious in some situations.
- **Observations:** Observations during assessment corroborated reports.
- **Diagnosis:** Utilizing this combined information to guide diagnosis, the clinician decided that Grant met criteria for a primary diagnosis of ADHD Predominantly Inattentive Presentation.



MHS
ASSESSMENTS



MHS
ASSESSMENTS

Progress Monitoring

Table 4.B. Conners CPT 3 Pre- and Post-Treatment Scores for Grant S.

Conners CPT 3 Scores	Pre-treatment Evaluation (Time 1)		Post-treatment Evaluation (Time 2)		Statistical Change
	F-score	Classification	F-score	Classification	
AD	62	Elevated	58	High Average	Time 1 = Time 2
Attention	74	Very Elevated	58	High Average	Time 1 > Time 2*
Commissions	72	Very Elevated	61	Elevated	Time 1 > Time 2*
Perseverations	65	Very Elevated	54	Average	Time 1 > Time 2*
WRT	68	Slow	58	A Little Slow	Time 1 > Time 2*
WRT SD	71	Very Elevated	61	Elevated	Time 1 > Time 2*
Variability	65	Elevated	60	Elevated	Time 1 = Time 2
WRT Block Change	52	Average	62	Elevated	Time 1 < Time 2
WRT ISI Change	52	Average	63	Elevated	Time 1 < Time 2

Note. The "*" and "**" symbols indicate scores that are statistically significant ($p < .10$) and/or at least 10 F-score points apart. Statistically significant changes are denoted by the "*" symbol.

MHS
ASSESSMENTS

CATA Report

- Very similar structure to the CPT 3 report
- Offers additional information about auditory laterality and mobility

CONNER'S
CPT3
CATA

MHS
ASSESSMENTS


K-CPT 2

- For age 4-7
- 7.5 minutes; 200 trials + 1 dummy trial
- Pictures of objects familiar to young children.
- 75% targets (everything except soccer ball)
- Presentation speed (Inter Stimulus Interval) can vary: 1.5 or 3.0 seconds
- Results can be broken down into blocks: 5 blocks with 40 trials each
- Dimensions of Attention Measured:
 - 1) Inattentiveness
 - 2) Impulsivity
 - 3) Sustained Attention (new)
 - 4) Vigilance (new)

CONNER'S
CPT3
CATA

MHS
ASSESSMENTS

Old K-CPT: Black Background



CONNER'S
CPT3
CATA

MHS
ASSESSMENTS

K-CPT 2: White background



CONNERA
CPT3
CATA

MHS
ASSESSMENTS

K-CPT 2 vs. CPT 3

	K-CPT 2	CPT 3
Admin Time	7.5 minutes	14 minutes
Stimuli	Pictures of common objects	letters
ISIs	1.5 & 3 seconds	1, 2, & 4 seconds
Stimuli Display Time	500ms	250ms

CONNERA
CPT3
CATA

MHS
ASSESSMENTS

Questions
about the
CPT3/CATA/
KCTP2?



References

- Conners, K. (2014). *Conners Continuous Performance Test 3rd Edition and Conners Continuous Auditory Test of Attention: Manual*. Toronto, Ontario, Canada: Multi-Health Systems.
- Dawson, P. (2014). Best practices in assessing and improving executive skills. In P. Harrison, & A. Thomas, *Best practices in school psychology student-level services* (pp. 269-286). Bethesda, MD: National Association of School Psychologists.
- Goldstein, S. & Naglieri, J. A. (2013). *Comprehensive Executive Function Inventory (CEFI): Technical manual*. Toronto, Ontario, Canada: Multi-Health Systems.
- Tobin, R., Scheider, W. & Landau, S.(2014). Best practices in the assessment of youth with attention deficit hyperactivity disorder within a multitiered services framework. In P. Harrison, & A. Thomas, *Best practices in school psychology data-based and collaborative decision making* (pp. 391-404). Bethesda, MD: National Association of School Psychologists.